BC 30

| | | |
|------------|--|--------------|
| D | Blutdruckmessgerät Gebrauchsanweisung | 2-11 |
| GB | Blood pressure monitor Instructions for use | 12-20 |
| F | Tensiomètre Mode d'emploi | 21-29 |
| E | Tensiómetro Manual de instrucciones | 30-38 |
| I | Misuratore di pressione Istruzioni per l'uso | 39-47 |
| TR | Bilgisayarlı tansiyon ölçer Kullanım kılavuzu | 48-56 |
| RUS | Прибор для измерения артериального давления на запястье Инструкция по применению | 57-67 |
| PL | Cisnieniomierz Instrukcja obsługi | 68-76 |
| | Electromagnetic Compatibility Information | 77-79 |

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Dear Customer,

thank you for choosing one of our products. Our name stands for high-quality, thoroughly tested products for the applications in the areas of heat, weight, blood pressure, body temperature, pulse, gentle therapy, massage and air. Please read these instructions for use carefully and keep them for later use, be sure to make them accessible to other users and observe the information they contain.

Best regards,
Your Beurer Team

1. Getting to know your instrument

The wrist blood pressure monitor is used for non-invasive measurement and monitoring of adults' arterial blood pressure. You can use it to measure your blood pressure quickly and easily, storing the results and displaying the progression of readings. A warning is issued for anyone suffering from cardiac arrhythmia.


The values determined are classified and graphically evaluated according to WHO guidelines. Keep these instructions carefully for further use and also let other users have access to them.

2. Important information

Signs and symbols

The following symbols are used in these instructions for use, on the packaging and on the type plate for the device and accessories:

| | |
|---|---|
|  | Caution |
|  | Note Note on important information |
|  | Follow instructions for use |
|  | Type BF applied part |
|  | Direct current |
|  | Disposal in accordance with EC Directive 2002/96/EC – WEEE (Waste Electrical and Electronic Equipment). |
|  | Manufacturer |

| | |
|---|--|
|  | Temperature limitation |
|  | Permissible storage humidity |
|  | Keep dry |
| SN | Serial number |
|  | The CE labelling certifies that the product complies with the essential requirements of Directive 93/42/EEC on medical products. |

Advice on use

- In order to ensure comparable values, always measure your blood pressure at the same time of day.
- Before every measurement, relax for about five minutes.
- If you want to perform several measurements on the same person, wait five minutes between each measurement.
- Do not take a measurement within 30 minutes after eating, drinking, smoking or exercising.
- Repeat the measurement if you are unsure of the measured value.
- The measurements taken by you are for your information only – they are not a substitute for a medical examination! Discuss the measurements with your doctor, and never base any medical decisions on them (e.g. medicines and their administration)!

- Do not use the blood pressure monitor on newborns, pregnant women or patients with preeclampsia.
- In the case of restricted circulation on the arm as a result of chronic or acute vascular diseases (including vascular constriction), the accuracy of the wrist measurement is limited. In this case you should avoid using an upper arm blood pressure monitor.
- Cardiovascular diseases may lead to incorrect measurements or have a detrimental effect on measurement accuracy. The same also applies to very low blood pressure, diabetes, circulatory disorders and arrhythmias as well as chills or shaking.
- The blood pressure monitor must not be used in connection with a high-frequency surgical unit.
- Only use the unit on people who have the specified wrist measurement for the device.
- Please note that when inflating, the functions of the limb in question may be impaired.
- During the blood pressure measurement, blood circulation must not be stopped for an unnecessarily long time. If the device malfunctions, remove the cuff from the arm.
- Do not allow sustained pressure in the cuff or frequent measurements. The resulting restriction of the blood flow may cause injury.
- Ensure that the cuff is not placed on an arm in which the arteries or veins are undergoing medical treatment, e.g. intravascular access or therapy, or an arteriovenous (AV) shunt.
- Do not use the cuff on people who have undergone a mastectomy.

- Do not place the cuff over wounds as this may cause further injury.
- The blood pressure monitor can only be operated with batteries.
- To conserve the batteries, the monitor switches off automatically if no buttons are pressed for one minute.
- The device is only intended for the purpose described in these instructions for use. The manufacturer is not liable for damage resulting from improper or careless use.

Storage and Care

- The blood pressure monitor is made up of precision electronic components. Accuracy of readings and the instrument's service life depend on careful handling.
 - You should protect the device from impact, moisture, dirt, major temperature fluctuations and direct exposure to the sun's rays.
 - Never drop the device.
 - Do not use near strong electromagnetic fields, i.e. keep it away from any radio systems and mobile phones.
- Do not press any buttons until the cuff is in position.
- If the instrument is not used for any length of time, we recommend removing the batteries.

Advice on batteries

- Batteries can be fatal if swallowed. You should therefore store the batteries and products where they are inaccessible to small children. If a battery has been swallowed, call a doctor immediately.

- Batteries should not be charged or reactivated with any other means, nor should they be taken apart, thrown in the fire or short-circuited.
- Remove the batteries from the instrument if they are worn out or if you are not going to use the instrument for any length of time. This prevents any damage as a result of leakage. Always replace all the batteries at the same time.
- Never use different types of battery, battery brands or batteries with different capacities. You should preferably use alkaline batteries.

Repair and disposal

- Batteries do not belong in domestic refuse. Used batteries should be disposed of at the collection points provided.
- Never open the instrument. If these instructions are not heeded, the warranty will be null and void.
- Never attempt to repair the instrument or adjust it yourself. We can no longer guarantee perfect functioning if you do.
- Repairs may only be performed by Customer Service or authorized dealers. However, always check the batteries and replace them if necessary prior to making any complaint.
- The appliance should be disposed of according to Regulation 2002/96/EC-WEEE (Waste Electrical and Electronic Equipment). In case of queries, please contact the municipal authorities responsible for waste disposal in your area.


3. Unit description

1. Display
2. Wrist cuff
3. WHO indicator
4. On/Off button 
5. Memory button **M**
6. Function button **SET**
7. Setting button **MODE**
8. Battery compartment lid


Information on the display:

1. Systolic pressure
2. WHO classification
3. Diastolic pressure
4. Calculated pulse value
5. Cardiac arrhythmia symbol 

6. Pulse symbol 
7. Memory space number
8. Battery replacement symbol 
9. Unit in mmHg
10. Memory display, average value
11. Time and date
12. User memory *no. 1, no. 2, no. 3*

4. Prepare measurement

Inserting battery

- Remove the battery cover (left side of the unit).
- Insert two 1.5V micro batteries (alkaline, type LR03). Make absolutely sure that you insert the batteries with the correct polarity as marked. Do not use rechargeable batteries.
- Replace the battery cover carefully.

If the “change battery” icon  is displayed, measurement is no longer possible and you must replace all the batteries. Spent batteries and rechargeable batteries do not constitute normal household waste! They are considered to be toxic waste and, as such, should be disposed of in special containers, toxic waste collection points or brought to electrical goods dealers.

Note: You will find these markings on batteries containing harmful substances: Pb = battery containing lead, Cd = battery containing cadmium, Hg = battery containing mercury.

Setting the date and time

This menu allows you to set the following functions, one after another:

Year → Month → Date → Time


It is essential to set the date and time. Otherwise, you will not be able to save your measured values correctly with a date and time and access them again later.

The time is displayed in the 24-hour format.

-  If you press and hold the function button **SET**, you can set the values more quickly.
- Press the setting button **MODE**. The year flashes on the display.
- Set the year with the function button **SET** and confirm with the setting button **MODE**.
- Set the month, day, hour and minutes and confirm each setting with the setting button **MODE**.
-  If you wish to change the date and/or time after it has already been set, press the setting button **MODE** twice in order to return to settings mode.

5. Measuring blood pressure

Positioning cuff

- Bare your left wrist, making sure that the circulation in the arm is not restricted by any clothes etc. that are too tight. Position the cuff on the inside of your wrist.
- Fasten the cuff with the Velcro fastening so that the upper edge of the monitor is positioned approx. 1 cm below the ball of your thumb.
- The cuff has to be fitted tightly around the wrist but should not constrict it.

Correct posture

- Rest for approx. 5 minutes before each measurement. Otherwise there may be divergences.
- You can perform the measurement while sitting or lying. To carry out a blood pressure measurement, make sure you are sitting comfortably with your arms and back leaning on something. Do not cross your legs. Place your feet flat on the ground. Make sure to rest your arm and move it. Always make sure that the cuff is at heart level. Otherwise significant deviations can occur. Relax your arm and the palm of your hand.

- In order not to distort the result, it is important to keep still during the measurement and not talk.

Performing the blood pressure measurement

- As described above, attach the cuff and adopt the posture in which you want to perform the measurement.
 - Start the blood pressure monitor with the on/off button ①. The most recently used user memory (n.a. 1, n.a.2 or n.a.3) appears. To change the user memory, press the function button **SET**.
 - Now start the measurement by holding the on/off button ①.
 - Following the display check, during which all numbers light up, the cuff inflates automatically. Whilst the cuff is pumping up, the device already calculates measurements for estimating the necessary pump pressure. If this pressure is insufficient, the device automatically pumps more pressure.
 - Then the pressure in the cuff is slowly released and as soon as a pulse is detected, the pulse symbol ♥ flashes.
- ① Measuring can be cancelled at any time by pressing the on/off button ① and the device switches to standby mode.
- Systolic pressure, diastolic pressure and pulse readings are displayed.
 - The measurement is automatically stored.
 - The symbol E_ appears if the measurement could not be performed properly. Observe the chapter “Error messages/ troubleshooting” in these instructions for use and repeat the measurement.
 - The device switches off automatically after 1 minute.

Wait at least 5 minutes before taking another measurement!

Evaluating results

Cardiac arrhythmia:

This instrument can identify possible cardiac arrhythmia disorders during measurement and if necessary indicates the measurement with the flashing icon ♥.

This may be an indicator for arrhythmia. Arrhythmia is a condition where the heart rhythm is abnormal as a result of defects in the bioelectrical system controlling the heart beat. The symptoms (omitted or premature heart beats, slow or excessively fast heart rate) may be caused, among other things, by heart disease, age, physical predisposition, excessive use of stimulants, stress or lack of sleep. Arrhythmia can only be ascertained through examination by your doctor. Repeat the measurement if the flashing icon ♥ is displayed after the measurement. Please note that you should rest for 5 minutes between measurements and not talk or move during the measurement. If the icon ♥ appears often, please contact your doctor. Any self-diagnosis and treatment based on the test results may be dangerous. It is vital to follow your doctor's instructions.

WHO classification:

According to WHO Guidelines/Definitions and the latest findings, the test results can be classified and evaluated according to the following chart:

| Range of blood pressure values | Systolic (in mmHg) | Diastolic (in mmHg) | Measure |
|--------------------------------|--------------------|---------------------|-------------------------------------|
| Grade 3: Severe hypertension | >=180 | >=110 | Seek medical advice |
| Grade 2: Moderate hypertension | 160-179 | 100-109 | Seek medical advice |
| Grade 3: Mild hypertension | 140-159 | 90-99 | Have it checked regularly by doctor |
| High-normal | 130-139 | 85-89 | Have it checked regularly by doctor |
| Normal | 120-129 | 80-84 | Check it yourself |
| Optimal | <120 | <80 | Check it yourself |

Source: WHO, 1999

The bar graph in the display and the scale on the unit indicate the range of the blood pressure which has been recorded. If the values for systolic and diastolic pressure are in two different WHO ranges (e.g. systolic in the “high-normal” range and diastolic pressure in the “normal” range) the graphic WHO classification on the unit indicates the higher range (“high-normal” in the example described).

6. Measurements

Save

- The results of every successful measurement are stored together with the date and time. There are 3 user memories each with 40 memory spaces. If there are more than 40 measurements, the oldest measurements are lost.

Display

- Select the desired user memory (m_0 , 1 , m_02 or m_03) by pressing the function button **SET**.
- Press the memory button **M**. The average value of the last three measurements is displayed first **AVG. 3**.

- If you continue to press the memory button **M**, the most recent individual measured values are displayed in turn with the date and time.

Delete

- To clear the memory of the relevant user memory, you must first select a user memory. Start a request for either average values or individual measured values and press the setting button **MODE** and the function button **SET** at the same time.
- All the values in the current user memory are deleted. $\square \square \square$ appears on the display.

The device switches off automatically after 1 minute.

7. Error message/troubleshooting

In the event of errors, the error message E_{-} appears on the display.

Error messages may appear in the following cases:

- EE : the cuff is incorrectly fastened, or you move or speak during the measurement,

- E : errors during measuring,
- E_2 : the pump pressure is higher than 300 mmHg,
- E_3 : errors in the recorded measurement data,
- EP : there is a system or unit error.

In such cases, repeat the measurement. Ensure that you do not move or speak.

If necessary, reinsert or replace the batteries.

8. Cleaning and Care

- Clean your blood pressure monitor carefully using a slightly damp cloth only.
- Do not use any detergents or solvents.
- Never hold the instrument under water as otherwise liquid can penetrate and damage the instrument.
- Never place any heavy objects on the instrument.

9. Technical details

| | |
|------------------------|---|
| Model no. | BC 30 |
| Measurement method | Oscillometric, non-invasive blood pressure measurement on the wrist |
| Measurement range | Cuff pressure 0–300 mmHg, systolic 40–280 mmHg, diastolic 40–280 mmHg, Pulse 40–199 beats/minute |
| Display accuracy | Systolic ± 3 mmHg, diastolic ± 3 mmHg, pulse $\pm 5\%$ of the value shown |
| Measurement inaccuracy | Max. permissible standard deviation according to clinical testing: systolic 8 mmHg/diastolic 8 mmHg |

| | |
|----------------------------------|--|
| Memory | 3 x 40 memory spaces |
| Dimensions | L 70 mm x W 72 mm x H 28,6 mm |
| Weight | Approx. 112 g (without batteries) |
| Cuff size | 135 to 195 mm |
| Permissible operating conditions | +5 °C to +40 °C, 15–93 % relative air humidity (non-condensing) |
| Permissible storage conditions | -25 °C to +70 °C, $\leq 93\%$ relative air humidity, 700–1060 hPa ambient pressure |
| Power supply | 2 x 1,5V  AAA batteries |
| Battery life | For approx. 250 measurements, depending on the blood pressure level and/or pump pressure |
| Accessories | Instruction for use, 2 x 1.5V AAA batteries, storage box |
| Classification | Internal supply, IPX0, no AP or APG, continuous operation, type BF applied part |

Technical information is subject to change without notification to allow for updates.

- This unit is in line with European Standard EN 60601-1-2 and is subject to particular precautions with regard to electromagnetic compatibility (EMC). Please note that portable and mobile HF communication systems may interfere with this unit. More details can be requested from the

stated Customer Service address or found at the end of the instructions for use.

- This device is in line with the EU Medical Devices Directive 93/42/EC, the “Medizinproduktegesetz” (German Medical Devices Act) and the standards EN 1060-1 (non-invasive sphygmomanometers, Part 1: General requirements), EN 1060-3 (non-invasive sphygmomanometers, Part 3: Supplementary requirements for electro-mechanical blood pressure measuring systems) and IEC 80601-2-30 (Medical electrical equipment – Part 2–30: Particular requirements for the safety and essential performance of automated non-invasive blood pressure monitors).
- The accuracy of this blood pressure monitor has been carefully checked and developed with regard to a long useful life. If using the device for commercial medical purposes, it must be regularly tested for accuracy by appropriate means. Precise instructions for checking accuracy may be requested from the service address.

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Electromagnetic Compatibility Information

Table 1

| <i>Guidance and manufacturer's declaration – electromagnetic emissions</i> | | |
|---|-------------------|--|
| The BC 30 is intended for use in the electromagnetic environment specified below. The customer or the user of the BC 30 should assure that it is used in such an environment. | | |
| Emissions test | Compliance | Electromagnetic environment – guidance |
| RF emissions CISPR 11 | Group 1 | The BC 30 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment. |
| RF emissions CISPR 11 | Class B | The BC 30 is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes. |
| Harmonic emissions IEC 61000-3-2 | Not applicable | |
| Voltage fluctuations/flicker emissions IEC 61000-3-3 | Not applicable | |

Table 2

| <i>Guidance and manufacturer's declaration – electromagnetic immunity</i> | | | |
|---|-----------------------------|----------------------------|--|
| The BC 30 is intended for use in the electromagnetic environment specified below. The customer or the user of the BC 30 should assure that it is used in such an environment. | | | |
| Immunity test | IEC 60601 test level | Compliance level | Electromagnetic environment – guidance |
| Electrostatic discharge (ESD) IEC 61000-4-2 | ±6 kV contact ±8 kV air | ±6 kV contact ±8 kV air | Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %. |
| Power frequency (50/60 Hz) magnetic field IEC 61000-4-8 | 3 A/m | 3 A/m | Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment. |

Table 3

| Guidance and manufacturer's declaration – electromagnetic immunity | | | |
|--|---|------------------|---|
| The BC 30 is intended for use in the electromagnetic environment specified below. The customer or the user of the BC 30 should assure that it is used in such an environment. | | | |
| Immunity test | IEC 60601 test level | Compliance level | Electromagnetic environment – guidance |
| Conducted RF IEC 61000-4-6 | 3V _{rms} 150 kHz to 80 MHz | Not applicable | <p>Portable and mobile RF communications equipment should be used no closer to any part of the BC 30, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance:</p> $d = 1.2 \sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ |
| Radiated RF IEC 61000-4-3 | 3V/m 80 MHz to 2.5 GHz | 3V/m | $d = 2.3 \sqrt{P} \quad 800 \text{ MHz to } 2.5 \text{ GHz}$ <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^a should be less than the compliance level in each frequency range.^b</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol: </p> |
| <p>NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.</p> <p>NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p> | | | |
| <p>^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the BC 30 is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the BC 30.</p> <p>^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.</p> | | | |

Table 4**Recommended separation distances between portable and mobile RF communications equipment and the BC 30**

The BC 30 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the BC 30 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the BC 30 as recommended below, according to the maximum output power of the communications equipment.

| Rated maximum output power of transmitter (W) | Separation distance according to frequency of transmitter (m) | | |
|---|---|---|--|
| | 150 kHz to 80 MHz Not applicable | 80 MHz to 800 MHz $d = 1.2 \sqrt{P}$ | 800 MHz to 2.5 GHz $d = 2.3 \sqrt{P}$ |
| 0.01 | – | 0.12 | 0.23 |
| 0.1 | – | 0.38 | 0.73 |
| 1 | – | 1.2 | 2.3 |
| 10 | – | 3.8 | 7.3 |
| 100 | – | 12 | 23 |

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.